

Residual Seal Force: A Powerful Vial Seal Quality Test



PDA
TRAINING

Carolina González Gaitán, PhD
Parenteral Packaging Scientist
Genesis Packaging Technologies

Agenda

1. Seal quality tests
2. Characterizing a “well-sealed” vial
3. Residual Seal Force – *Concept, basis of testing, methodology, variability considerations, significance and use of RSF test*
4. Studies – *Effect of time, effect of FO button, correlation with CCIT*
5. Takeaways

Seal Quality Tests

- USP <1207.3>:


*“Package seal quality tests” are checks used to **characterize and monitor** the **quality and consistency** of a seal parameter related to the package seal, providing some assurance of the package’s ability to remain integral*

Quality tests ensure that **seal attributes, package materials, package components** and/or the **assembly process** are consistently kept within established limits, thus further supporting **package integrity**

Seal quality tests are **NOT** leak tests


“Well-Sealed” Vial

- Sufficient compression to achieve Leak Rate Cut-off
- An applied force compresses the stopper flange.
 1. Cross section of the component(s)
 2. Durometer (hardness) of the rubber
 3. The percent of compression required to achieve leak rate cut-off



“Well-Sealed” Vial

- Sufficient compression to achieve Leak Rate Cut-off
- An applied force compresses the stopper flange.
 1. Cross section of the component(s)
 2. Durometer (hardness) of the rubber
 3. The percent of compression required to achieve leak rate cut-off




Morton, Dana K. "Quantitative and Mechanistic Measurements of Parenteral Vial Container/Closure Integrity. Leakage Quantitation" *PDA J of Pharm Sci and Technol* 1989, 43 (2) 88-97

Residual Seal Force - RSF

- RSF is the strain a compressed elastomeric rubber stopper flange continues to exert on the vial sealing surface after the crimping of an aluminum seal
- RSF is an easy-to-use quantitative method to standardize seal quality regardless of the capping equipment used for crimping
- RSF helps to set up capping parameters to ensure consistency and ease capper validations
- Correlation of RSF with CCITs will provide guidance on setting acceptable ranges



Basis of RSF Testing




- Upon capping, the stopper flange is compressed against the vial land sealing surface
- The stopper flange acts like a “compressed spring”
- The tester apply a force on the cap and stopper
- When the tester force exceeds the closure compression force → RSF

R. Mathaes et al. “The pharmaceutical vial capping process: Container closure systems, capping equipment, regulatory framework, and seal quality tests” *European Journal of Pharmaceutics and Biopharmaceutics* 99 (2016) 54–64

RSF Tester and Methodology



Determining RSF



- Stress-strain curve (green) is a combination of the viscous and elastic response to the stress from tester load
- RSF is determined using the stress-strain curve: the “knee” (yellow)
- An algorithm* is applied, using the 1st (purple) and 2nd (blue) derivatives to accurately identify that knee

* Ludwig J, Nolan P, Davis C, Automated method for determining Instron residual seal force of glass vial/rubber stopper closure systems, *PDA J of Pharm Sci and Technol* 1993, 47 (5) 211-253

Variability Considerations



Instrument

- Gage R&R
- Custom compressed spring fixture ~2%
- Orientation & centering
- Anvil design
- Button removal



Components Variation

- Dimensional tolerances
- Stack-up, interference fit
- Mismatch of components



Time

- Elastomer relaxation



Capping Process

- Optimization of settings
- One head vs. Multiple heads

Significance & Use of RSF Method



Package Development

- Determine effects of CCS component variables
- Characterize a “well-sealed” vial



Validation

- Establish optimum capping parameters
- Evaluate variation



Production


- Verify capping equipment set-up
- Capping process monitor



RSF – Influence of Time

Influence of Elastomer Relaxation

Maxwell-Weichert
Degeneration Curve




- Elastomer is the base material of the stopper
- Exhibit viscoelastic behavior
- Relaxes over time → RSF decay over time

Morton D., Lordi N. "Residual Seal Force Measurements of Parenteral Vials: I. Methodology" PDA J Pharm Sci and Technol 1998, 42 23-29

RSF – Time


Figure 1: RSF and helium leak testing data for vial CCS using a 20 mm butyl elastomer stopper and a 10 mL glass vial fully filled with helium at ambient pressure, tested at ambient temperature through a vacuum chamber [8, 9].



Zeng, Q. “Critical Time- & Temperature- Dependent Container Closure Integrity (CCI) Through the Sealed Drug Product Life Cycle” PDA Parenteral Packaging Conference, Rome, Italy; 2018

RSF – Time

Time dependent RSF testing at ambient conditions and modeling fit
(20mm serum stopper, seal, & vial)



Zeng, Q.; Zhao, C.; "Time-Dependent Testing Evaluation and Modeling for Rubber Stopper Seal Performance." *PDA J Pharm Sci and Tech* 2018, 72 134-148

RSF – Time

Statistical Data Generated of 20 Vials from the RSF Time Course

Time	Mean RSF (N) (n = 20)	Difference in Mean	RSD%
1 minute	62.7	-	9.9
10 minutes	54.0	8.7	11.0
90 minutes	53.1	0.9	7.0
1 day	52.1	1.0	9.6
7 days	51.0	0.9	11.1
21 days	50.5	0.5	10.2

Adapted from: Ovadia, R; Streubel, A; et al. “Quantifying the Vial Capping Process: Residual Seal Force and Container Closure Integrity” *PDA J of Phar Sci and Technol*, 2019 73 (1) 2-15


RSF – Time

- Stress-relaxation of the rubber stopper is time-dependent affecting the sealing force
- Rubber will relax with time
 - RSF decay
 - Greater variability at $t < 10$ min
 - Greater decrease with higher crimping forces


RSF - Flip-Off Cap Impact

Flip-Off Cap Impact


Without flip-off button




One clear minimum on 2nd derivative



With flip-off button




More complex and noisier signal



R. Mathaes et al. "Impact of Vial Capping on Residual Seal force and Container Closure Integrity" *PDA J Pharm Sci and Tech* 2016, 70 12-29


Flip-Off Cap Impact

Without flip-off button



Low variability
Distinctive RSF groups

With flip-off button



High variability
Difficult to distinguish among RSF groups

R. Mathaes et al. "Impact of Vial Capping on Residual Seal force and Container Closure Integrity" *PDA J Pharm Sci and Tech* 2016, 70 12-29


RSF – FO Cap

- The flip-off button adds complexity to the system, preventing a clean transition of the force applied by the RSF tester
 - The stress-strain curve is more complex – sometimes with 2 minima
 - Higher variability
- More reliable results without the flip-off button → Destructive

Correlation with CCIT

Correlation - RSF to Compression


- CCS:
 - 10R Vial
 - 20 mm serum stopper
- Sealing parameters:
 - Four (4) crimping pressures / RSF targets
- Compression, RSF and He leak




Correlation - RSF to He Leak Rate

- Kirsch criterion*: Helium leak rates lower than 6×10^{-6} std cc/s have been associated with acceptable microbial challenge results
- Low group have several samples that failed based on the Kirsch Criterion

*Kirsch, L et al. "Pharmaceutical container/closure integrity II: The relationship between microbial ingress and helium leak rates in rubber-stoppered glass vials" *PDA J of Pharm Sci and Technol* 51 (5) 195-202 (1997)




Correlation - RSF to HVLD



S. Orosz and D Guazzo, "Leak Detection and Product Risk Assessment" presented at PDA Annual Meeting, Mar 2010, Orlando, FL


Correlation - RSF to HVLD



S. Orosz and D Guazzo, "Leak Detection and Product Risk Assessment" presented at PDA Annual Meeting, Mar 2010, Orlando, FL

Correlation - RSF to HSA

- CCS:
 - 2 ml Vial EU BB, 13 mm Serum Stopper
 - Five (5) vial stopper combinations (A – E)
- Sealing parameters:
 - Three (3) crimping pressures – RSF targets
- Storage:
 - Four (4) storage temperatures




Duncan, D.; Asselta, R. "Correlating Vial Seal Tightness to Container Closure Integrity at Various Storage Temperatures" proceedings of PDA Parenteral Packaging Conference, Frankfurt, Germany; (2015)

Correlation - RSF to HSA

At -80°C:

- *Package A*: 24% failures at low compression setting
- *Package B*: 7% failures at low compression setting
- *Package C*: 0% failures at low compression setting, 4% failures at Nominal compression setting
- *Package D*: 10% failures at low compression setting
- *Package E*: 4% failures at low compression setting



Duncan, D.; Asselta, R. "Correlating Vial Seal Tightness to Container Closure Integrity at Various Storage Temperatures" proceedings of PDA Parenteral Packaging Conference, Frankfurt, Germany; (2015)

RSF – CCIT

- Correlation of RSF to CCITs will provide guidance on setting acceptable ranges
- Once optimal RSF range is established, it can be used to standardize seal quality regardless the capping equipment used for crimping

Takeaways

- RSF is a reliable and precise measurement to assess the quality of sealed vial and predict CCI failure
- The stopper compression is a function of RSF
- Correlation of RSF and CCITs provides guidance on setting acceptable ranges, allowing comparison among different capping equipment & sites



Genesis Packaging Technologies

CONNECTING
PEOPLE
SCIENCE AND
REGULATION®



Website:

www.gen-techno.com



Email:

info@gen-techno.com

CGonzalez@gen-techno.com



Phone:

+1 800 552 9980

+1 613 294 9203

Thank you!



PDA
TRAINING

pda.org