



Process Validation for Lyophilized Drug Products: Comparing a Program for Continued Process Verification in Different Lyophilized Products

Abirami Natesh and Denise Miller

Lyophilization Technology, Inc.

INTRODUCTION

Maintaining process control in lyophilization is essential for ensuring product quality and regulatory compliance. The FDA’s 2011 guidance on process validation outlines a three-stage approach to ensuring a state of control throughout the product lifecycle: Stage 1 (Process Design) establishes critical parameters and proven acceptable ranges (PARs), Stage 2 (Process Qualification) validates these conditions, and Stage 3 (Continued Process Verification, CPV) ensures ongoing monitoring and control.

While Critical Quality Attributes (CQAs) can be effectively trended and reported annually, monitoring Critical Process Parameters (CPPs)—such as shelf temperature and chamber pressure—presents greater challenges, particularly in multi-product facilities. Time, however, is typically excluded as a monitored parameter since freeze-dryers operate with pre-programmed time settings that are qualified with the unit itself, eliminating the need for ongoing trending. Traditional statistical methods for process variability analysis, though effective for single-product evaluations, become inefficient when applied across multiple lyophilized products.

A standardized approach using variation from setpoint as a process control metric offers a more efficient and scalable solution. By assessing deviations relative to predefined targets or PARs, manufacturers can streamline data analysis, enhance process consistency, and simplify ongoing monitoring without relying on product-specific statistical models. This study evaluates the effectiveness of this methodology across multiple lyophilized products, demonstrating its potential to improve CPV, optimize operational efficiency, and ensure regulatory compliance in multi-product lyophilization environments.

VARIATION FROM SETPOINT AS A PROCESS CONTROL METRIC

Tracking deviations from setpoints in parameters like shelf temperature and chamber pressure offers significant insights for process monitoring by identifying trends that could impact product quality, allowing for proactive corrective actions. This approach improves comparability and efficiency across products and freeze-drying cycles by providing a consistent metric for evaluating process performance. Since variation from setpoint is independent of batch cycles and product types, it is adaptable and scalable, ensuring effective continued process verification. In multi-product facilities with varying lyophilization cycles, using this metric standardizes data analysis, while control charts enhance process verification by visually assessing stability and detecting deviations.

DATA

The data in Figures 1-5 was compiled from two different products, Product A and B each having different lyophilization cycles. Throughout the nine batches, both cycles had additional refinements which modified their parameters and product load. However, since this method tracks variation from setpoint and not specific values it can be used for any cycle even as the parameters change. In the graphs, the minimum, maximum, and average variation for each batch are shown.

These graphs show that during freezing, primary drying, and secondary drying, the cycles were within the allowable variation within the PAR from the shelf temperature setpoint of $\pm 5^{\circ}\text{C}$.




Figure 1: Freezing Step Minimum, Maximum, and Average Temperature Variation across multiple lots and multiple products




Figure 2: Primary Drying Step Minimum, Maximum, and Average Temperature Variation across multiple lots and multiple products




Figure 3: Secondary Drying Step Minimum, Maximum, and Average Temperature Variation across multiple lots and multiple products

DATA (continued)

Figures 1-3 show how the variation from setpoint allows for ease of reviewing data for the multiple batches at each individual step of the process. To further compress the data to allow for ease of trending, the entire cycle can be compiled using variation from the setpoint as seen in Figure 4.




Figure 4: Minimum, Maximum, and Average Temperature Variation across the Lyophilization Cycle

For chamber pressure, the allowable variation range is greater than that of temperature; in this case, the range of variation is $\pm 20\mu\text{m}$. It is logical to chart chamber pressure data separately to easily trend the data. Although chamber pressure is recorded throughout the entire cycle, it is only during primary and secondary drying that the chamber pressure becomes a critical process parameter. This is because vacuum is typically pulled during these stages and thus requires monitoring to ensure chamber pressure stays within the PAR.




Figure 5: Minimum and Maximum Chamber Pressure Variation during Primary and Secondary Drying compared to the Overall Average Variation

Using the chamber pressure control chart, it is simple to visually identify pressure excursions that occurred as seen in Figure 5.

CONCLUSION

Effective process control in lyophilization is essential for ensuring product quality, regulatory compliance, and operational efficiency, particularly in multi-product facilities. Traditional statistical methods, while effective for single-product evaluations, become impractical when applied across multiple lyophilized products with varying process parameters. By utilizing variation from setpoint as a process control metric, manufacturers can establish a standardized, scalable approach to using a control chart for the Critical Process Parameters (CPPs) of shelf temperature and chamber pressure. This method enhances process consistency, allows for proactive adjustments, and simplifies continued process verification (CPV) without the need for product-specific statistical models.

Integrating control charts further strengthens process monitoring by providing a visual representation of process stability across multiple batches and products. The case study findings confirm that this combined approach effectively tracks deviations, ensures processes remain within acceptable PAR boundaries, and facilitates informed decision-making. By adopting variation from setpoint and control charts as key process monitoring tools, manufacturers can improve efficiency, enhance product consistency, and maintain robust control over lyophilization cycles, ultimately supporting long-term process validation and regulatory compliance.

REFERENCES

Center for Drug Evaluation and Research. “Process Validation: General Principles and Practices.” U.S. Food and Drug Administration, FDA, www.fda.gov/regulatory-information/search-fda-guidance-documents/process-validation-general-principles-and-practices.

Bossert, K., Trappler, E., Sitapara, A. (2015). Process Validation for Lyophilized Drug Products: Developing a Program for Continued Process. Lyophilization Technology, Inc., Ivyland, PA

ACKNOWLEDGEMENTS

Thank you to the dedicated staff at Lyophilization Technology, Inc.

